

0057849

H1797

Lionville Laboratory, Inc.
INORGANIC ANALYTICAL DATA PACKAGE FOR
TNUHANFORD B02-059 H1797

DATE RECEIVED: 05/30/02

LVL LOT # :0205L803

CLIENT ID /ANALYSIS	LVL #	MTX	PREP #	COLLECTION	EXTR/PREP	ANALYSIS
---------------------	-------	-----	--------	------------	-----------	----------

B14R98

LEAD, TOTAL	001	W	02L0290	05/23/02	06/04/02	06/06/02
LEAD, TOTAL	001 REP	W	02L0290	05/23/02	06/04/02	06/06/02
LEAD, TOTAL	001 MS	W	02L0290	05/23/02	06/04/02	06/06/02

LAB QC:

LEAD LABORATORY	LC1 BS	W	02L0290	N/A	06/04/02	06/05/02
LEAD, TOTAL	MB1	W	02L0290	N/A	06/04/02	06/05/02



RECEIVED
AUG 12 2002
EDMC



Analytical Report

Client: TNU-HANFORD B02-059
LVL#: 0205L803
SDG/SAF#: H1797/B02-059


W.O.#: 11343-606-001-9999-00
Date Received: 05-30-02

METALS CASE NARRATIVE

1. This narrative covers the analysis of 1 water sample.
2. The sample was prepared and analyzed in accordance with methods checked on the attached glossary.
3. All analyses were performed within the required holding times.
4. Please refer to the Sample Receipt Check List for sample discrepancies in LvLI's sample acceptance policy.
5. All Initial and Continuing Calibration Verifications (ICV/CCVs) were within the 90-110% control limits.
6. All Initial and Continuing Calibration Blanks (ICB/CCBs) were within control limits (less than the PQL).
7. The preparation/method blank (MB) was within method criteria {less than the Practical Quantitation Limit (3X the IDL), MB value less than 5% of the RCRA limit, or samples greater than 20X MB value}. Refer to the Inorganics Method Blank Data Summary.
8. All ICP Interference Check Standards were within control limits.
9. The laboratory control sample (LCS) was within the 80-120% control limits. Refer to the Inorganics Laboratory Control Standards Report.
10. The matrix spike (MS) recovery was within the 75-125% control limits. Refer to the Inorganics Accuracy Report.
11. The duplicate analysis was within the 20% Relative Percent Difference (RPD) control limits. Refer to the Inorganics Precision Report.

The results presented in this report relate only to the analytical testing and conditions of the samples at receipt and during storage. All pages of this report are integral parts of the analytical data. Therefore, this report should only be reproduced in its entirety of 14 pages.

12. For the purposes of this report, the data has been reported to the Instrument Detection Limit (IDL). Values between the IDL and the Practical Quantitation Limit (PQL) are acquired in a region of less-certain quantification.
13. I certify that this sample data package is in compliance with SOW requirements, both technically and for completeness, other than the conditions detailed above. Release of the data contained in this hard-copy data package has been authorized by the Laboratory Manager or a designee, as verified by the following signature.


Iain Daniels
Laboratory Manager
Lionville Laboratory Incorporated
gmb/m05-803


Date

METALS METHOD GLOSSARY

The following methods are used as reference for the digestion and analysis of samples contained within this

Lot#: 02052803

Leaching Procedure: ☐ 1310 ☐ 1311 ☐ 1312 ☐ Other: _____

CLP Metals ☐ Digestion and Analysis Methods: ☐ ILM03.0 ☐ ILM04.0

Metals Digestion Methods: ☒ 3005A ☐ 3010A ☐ 3015 ☐ 3020A ☐ 3050B ☐ 3051 ☐ 200.7 ☐ SS17
☐ Other: _____

Metals Analysis Methods

	SW846	EPA	STD MTD	EPA OSWR	USATHAMA
Aluminum	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Antimony	<input type="checkbox"/> 6010B <input type="checkbox"/> 7041 ¹	<input type="checkbox"/> 200.7 <input type="checkbox"/> 204.2			<input type="checkbox"/> 99
Arsenic	<input type="checkbox"/> 6010B <input type="checkbox"/> 7060A ¹	<input type="checkbox"/> 200.7 <input type="checkbox"/> 206.2	<input type="checkbox"/> 3113B		<input type="checkbox"/> 99
Barium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Beryllium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Bismuth	<input type="checkbox"/> 6010B ¹	<input type="checkbox"/> 200.7 ¹		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Boron	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Cadmium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7131A ¹	<input type="checkbox"/> 200.7 <input type="checkbox"/> 213.2			<input type="checkbox"/> 99
Calcium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Chromium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7191 ¹	<input type="checkbox"/> 200.7 <input type="checkbox"/> 218.2			<input type="checkbox"/> SS17
Cobalt	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Copper	<input type="checkbox"/> 6010B <input type="checkbox"/> 7211 ¹	<input type="checkbox"/> 200.7 <input type="checkbox"/> 220.2			<input type="checkbox"/> 99
Iron	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Lead	<input checked="" type="checkbox"/> 6010B <input type="checkbox"/> 7421 ¹	<input type="checkbox"/> 200.7 <input type="checkbox"/> 239.2	<input type="checkbox"/> 3113B		<input type="checkbox"/> 99
Lithium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7430 ¹	<input type="checkbox"/> 200.7		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Magnesium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Manganese	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Mercury	<input type="checkbox"/> 7470A ¹ <input type="checkbox"/> 7471A ¹	<input type="checkbox"/> 245.1 ² <input type="checkbox"/> 245.5 ¹			<input type="checkbox"/> 99
Molybdenum	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Nickel	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Potassium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7610 ¹	<input type="checkbox"/> 200.7 <input type="checkbox"/> 258.1 ¹			<input type="checkbox"/> 99
Rare Earths	<input type="checkbox"/> 6010B ¹	<input type="checkbox"/> 200.7 ¹		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Selenium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7740 ¹	<input type="checkbox"/> 200.7 <input type="checkbox"/> 270.2	<input type="checkbox"/> 3113B		<input type="checkbox"/> 99
Silicon	<input type="checkbox"/> 6010B ¹	<input type="checkbox"/> 200.7		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Silica	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Silver	<input type="checkbox"/> 6010B <input type="checkbox"/> 7761 ¹	<input type="checkbox"/> 200.7 <input type="checkbox"/> 272.2			<input type="checkbox"/> 99
Sodium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7770 ¹	<input type="checkbox"/> 200.7 <input type="checkbox"/> 273.1 ¹			<input type="checkbox"/> 99
Strontium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Thallium	<input type="checkbox"/> 6010B <input type="checkbox"/> 7841 ¹	<input type="checkbox"/> 200.7 <input type="checkbox"/> 279.2 <input type="checkbox"/> 200.9			<input type="checkbox"/> 99
Tin	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Titanium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Uranium	<input type="checkbox"/> 6010B ¹	<input type="checkbox"/> 200.7 ¹		<input type="checkbox"/> 1620	<input type="checkbox"/> 99
Vanadium	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Zinc	<input type="checkbox"/> 6010B	<input type="checkbox"/> 200.7			<input type="checkbox"/> 99
Zirconium	<input type="checkbox"/> 6010B ¹	<input type="checkbox"/> 200.7 ¹		<input type="checkbox"/> 1620	<input type="checkbox"/> 99

Other: _____

Method: _____

METHOD REFERENCES AND DATA QUALIFIERS

DATA QUALIFIERS

U = Indicates that the parameter was not detected at or above the reported limit. The associated numerical value is the sample detection limit.

* = Indicates that the original sample result is greater than 4x the spike amount added.

ABBREVIATIONS

MB = Method or Preparation Blank.
MS = Matrix Spike.
MSD = Matrix Spike Duplicate.
REP = Sample Replicate
LCS = Laboratory Control Sample.
NC = Not calculated.

ANALYTICAL METAL METHODS

1. Not included in the method element list.
2. Modified Hg: Hg1 and Hg2 require less total volume of digestate due to the autosampler analysis. Sample volumes and reagents for mercury determinations in water and soil have been proportionately scaled down to adapt to this semi-automated technique. The sample volume used for water analysis is 33 mL. For soils, 0.1 grams of sample is taken to a final volume of 50 mL (including all reagents).
3. Modified Hg: Hg1 and Hg2 require less total volume of digestate due to the autosampler analysis. Sample volumes and reagents for mercury determinations in water and soil have been proportionately scaled down to adapt to this semi-automated technique. The sample volume used for water analysis is 33 mL. For soils, three 0.1 gram of sample is taken to a final volume of 50 mL (including all reagents).
4. Flame AA.
5. Graphite Furnace AA.

Lionville Laboratory, Inc.

INORGANICS DATA SUMMARY REPORT 06/13/02

CLIENT: TNUHANFORD B02-059 H1797

LVL LOT #: 0205L803

WORK ORDER: 11343-606-001-9999-00

SAMPLE	SITE ID	ANALYTE	RESULT	UNITS	REPORTING LIMIT	DILUTION FACTOR
-----	-----	-----	-----	-----	-----	-----
-001	B14R98	Lead, Total	1300	UG/L	1.7	1.0

Lionville Laboratory, Inc.

INORGANICS METHOD BLANK DATA SUMMARY PAGE 06/13/02

CLIENT: TNUHANFORD B02-059 H1797

LVL LOT #: 0205L803

WORK ORDER: 11343-606-001-9999-00

SAMPLE	SITE ID	ANALYTE	RESULT	UNITS	REPORTING LIMIT	DILUTION FACTOR
BLANK1	02L0290-MB1	Lead, Total	9.0	UG/L	1.7	1.0

Lionville Laboratory, Inc.

INORGANICS ACCURACY REPORT 06/13/02

CLIENT: TNUHANFORD B02-059 H1797

LVL LOT #: 0205L803

WORK ORDER: 11343-606-001-9999-00

SAMPLE	SITE ID	ANALYTE	SPIKED SAMPLE	INITIAL RESULT	SPIKED AMOUNT	%RECOV	DILUTION FACTOR (SPK)
-001	B14R98	Lead, Total	1890	1300	500	118.2	1.0

Lionville Laboratory, Inc.

INORGANICS PRECISION REPORT 06/13/02

CLIENT: TNUHANFORD B02-059 H1797

LVL LOT #: 0205L803

WORK ORDER: 11343-606-001-9999-00

SAMPLE	SITE ID	ANALYTE	INITIAL RESULT	REPLICATE RPD	DILUTION FACTOR (REP)
-----	-----	-----	-----	-----	-----
-001REP	B14R98	Lead, Total	1300	1310 0.89	1.0

Lionville Laboratory, Inc.

INORGANICS LABORATORY CONTROL STANDARDS REPORT 06/13/02

CLIENT: TNUHANFORD B02-059 H1797
 WORK ORDER: 11343-606-001-9999-00

LVL LOT #: 0205L803

SAMPLE	SITE ID	ANALYTE	SPIKED SAMPLE	SPIKED AMOUNT	UNITS	%RECOV
-----	-----	-----	-----	-----	-----	-----
LCS1	02L0290-LC1	Lead, LCS	2500	2500	UG/L	99.8

22056803

FIELD PERSONNEL: COMPLETE ONLY SHADED AREAS

Discrepancies Between
Samples Labels and
COC Record Y or N
NOTES:
011511013993

SAMPLE RECEIPT CHECKLIST

CLIENT: T/11U-Hanford

Purchase Order/Project:

DATE: 5.30-02

SAF# / SOW# / Release #: 802-059

Laboratory SDG #: 0205 L803

NOTE: ALL ENTRIES MARKED "NO" MUST BE EXPLAINED IN THE COMMENT SECTION

- | | | | | |
|--|---|--|---|---|
| 1. Custody seals on coolers or shipping container intact, signed and dated? | <input checked="" type="checkbox"/> Yes | <input type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |
| 2. Outside of coolers or shipping containers are free from damage? | <input checked="" type="checkbox"/> Yes | <input type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |
| 3. Airbill # recorded? | <input checked="" type="checkbox"/> Yes | <input type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |
| 4. All expected paperwork received (coc and other client specific: historical data, alpha/beta or other screening data as applicable)? (paperwork sealed in plastic bag and taped to inside lid) | <input checked="" type="checkbox"/> Yes | <input type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |
| 5. Sample containers are intact? | <input checked="" type="checkbox"/> Yes | <input type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |
| 6. Custody seals on sample containers intact, signed and dated? | <input checked="" type="checkbox"/> Yes | <input type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |
| 7. All samples on coc received? | <input checked="" type="checkbox"/> Yes | <input type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |
| 8. All sample label information matches coc? | <input type="checkbox"/> Yes | <input checked="" type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # <u>(1)</u> |
| 9. Laboratory QC samples designated on coc? (QC stickers placed on bottles?) | <input type="checkbox"/> Yes | <input type="checkbox"/> No | <input checked="" type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |
| 10. Shipment meets LVLJ Sample Acceptance Policy? (identify all bottles not within policy. See reverse side for policy) | <input type="checkbox"/> Yes | <input checked="" type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # <u>(1)</u> |
| 11. Where applicable, bar code labels are affixed to coc? | <input type="checkbox"/> Yes | <input type="checkbox"/> No | <input checked="" type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |
| 12. coc signed and dated? | <input checked="" type="checkbox"/> Yes | <input type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |
| 13. coc faxed or emailed to client? | <input type="checkbox"/> Yes | <input type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |
| 14. Project Manager/Client contacted concerning discrepancies? (name/date) | <input type="checkbox"/> Yes | <input type="checkbox"/> No | <input type="checkbox"/> N/A | <input type="checkbox"/> see Comment # |

Cooler # / temp and Comments:

SPK
1511013993

Laboratory Sample Custodian:

Carlo Huang

Laboratory Project Manager:

SAMPLE DIGESTION RECORD

Digestion Batch #: 0210290
 Date/Time Initiated: 06/04/02
 Date/Time Completed: 06/04/02
 Analyst(s): EMP
 Matrix: Soil Water Other: _____
 Instr. Type: AA (ICP)
 Parameters: SEE BACKLOG

Method: SW
 (circle)

SOP: L-SPI-3020 Rev. 00

3005A
 3010A
 3015
 3020A
 7060A (As/Se)
 7760A (Ag)
 3050B
 3051
 DW 200.7 (1994)
 200.9
 3113B
 MCAWW 200.7 (1982)
 200 (AA)
 206.2 (As/Se)
 SM 3030C (NC)

(Digested) Undigested (circle one)

Balance #: NIA

Balance Cal Verif: Y (NIA)

Hot Plate Temp: 92°C

CLP ILM03.0
 ILM04.0

Other _____

TN4/WESKEM

COC Batch #	Spike Vol(s) (mL)	Initial Wt/Vol (g/mL)	Final Vol (mL)	pH <	Type: To/So/ TC	Texture	Color/Appearance	Artifact	Turb
0205L733-001		100	100	42	Tot	NIA	Cloudy	NIA	NIA
-001R									
-001S	1.0								
0205L803-001							Clear		
-001R									
-001S	1.0								
0205L778-013									
0205L780-033									
-094									
0205L624-001									
0205L724-001									
0205L725-001									
0210290-MB1				NIA					
-LC1	1.0								
<div style="text-align: center;"> <u>Prep</u> <u>06-04-02</u> </div>									

Spiking IDs:

MS #: 8100-02-05
-06
-07
6072-48-14
 LCS #: 6072-49-05
-06
-07
-08

Reagent IDs:

HNO₃
 HCL
 H₂O₂
 1:1 HNO₃
 1:1 HCL
X08043
2612V27006
2004/02
8728-025-02
8727-025-03

File ID#: IL029001

LIMS Transfer: (Y) N updated

Data Review By/Date:

upp 6/7/02

14